www.ifu.oasi.care

EN Instruction for use

ASISU®

Perineal Protection Device

 oasicare ApS
 www.oasi.care
 support@oasi.care

Single-use only

ASISU® is a medical device intended for use by qualified birth attendants during the second stage of vaginal childbirth, particularly during the crowning phase. The device consists of an internal part positioned in the vaginal canal under the fetal head and an external part placed on the perineal skin, connected by a flexible neck at the posterior fourchette. The device provides mechanical support to the perineum and helps distribute forces during fetal head expulsion with the aim of reducing the risk of perineal tears. The device is supplied ready to use in a pouch allowing visual inspection before opening.

Intended Purpose (Extended statement):

This medical device is intended to protect the perineum during the second stage of vaginal childbirth. It is intended to be used by qualified midwives, obstetric clinicians and other professional birth attendants, in clinical and community-birth environments, to reduce the risk of perineal tears and associated adverse consequences.

Indications for use:

Use during the second stage of labor, specifically in the crowning phase.

Contraindications:

- The device is not intended for patients with known perineal abnormalities or complications requiring surgical intervention during delivery.
- The device is not intended for use in births requiring instrumental intervention.
- No side-effects observed upon the use of the device as intended.

Intended user:

- Use only by qualified midwives, obstetricians, or other trained healthcare professionals attending childbirth.
- Claimed clinical benefits: Reduce the risk of perineal tears and associated adverse consequences.
- Performance Characteristics: Flexible internal and external adhesion elements provide perineal support and reduce tissue stress peak.
- Patient target group: Pregnant individuals undergoing vaginal childbirth.

Requirement of:

1. Special facilities:

- Clinical birth environments that include hospital wards.
- Community birth environments that include birth centers and planned home births.

2. Training Video:

For additional guidance on device preparation and application, click [here](#) to view the training video.

⚠ Warning & Caution:

- The device should only be used once to avoid device failure or sub-optimal performance, and infections.
- Dispose of the device as waste per local regulations.

⚠ Precaution:

- Position the outer part according to the woman's anatomy so that it does not cover or contact the anus. If necessary, trim the device before application to ensure correct fit.
- Where episiotomy is clinically indicated, it takes priority over use of the device. The midwife shall ensure that the device does not interfere with performance of the episiotomy.

eIFU Version: Issue 3


Date of issue: 2026/04/07

If a serious incident occurs in relation to ASISU®, please report it without delay to oasicare ApS (email: support@oasicare) and to the competent authority of the EU Member State where you are established.

The device is used in the final part of the second stage of labor

The second stage of labor starts: Prepare device

Usage time: 10–20 min


1. Insert the internal part at the end of a contraction.
2. Gently guide it into position while the fetal head recedes during the pause between contractions.
3. Apply gel if needed to facilitate insertion.
4. As the fetal head advances, it may push the internal part outward. Maintain the position of the internal part by applying gentle manual counter-pressure.
5. Dry the perineum with a clean cloth.
6. Gently remove the protective liner from the outer part.
7. Place the outer part on the perineum.
8. Position the wings according to the woman's anatomy. Ensure that the device does not cover the anus.
9. Press the outer part firmly in place for 30 seconds to activate the adhesive.
10. Once the device is secured, standard techniques for perineal protection may be used.
11. When the fetal head is delivered, maintain the internal part in position with gentle manual pressure.

12. After delivery, gently remove the adhesive outer part from the skin.

13. Carefully withdraw the internal part.

14. Dispose of the device as regular waste.

Copyright © 2026 Oasicare ApS. All rights reserved.

This document is proprietary to Oasicare ApS. No part of this document may be reproduced, distributed, or transmitted in any form or by any means without prior written permission from Oasicare ApS.

European Conformity

Catalogue number

Medical device

Consult electronic instructions for use

Unique device identifier

This package contains 1 unit

Do not use if package is damaged and consult instruction for use

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Document: OC-000679-RC**Approved version: 3****Title: eIFU****Approved By:**

CTO - Alireza Kouchaki Tue Apr 7 12:08:48 GMT 2026

Approved & Signed 1IK0Do4AFr16aqx04YY5SOSXzFs

QA/RA - Abeer Prakash Singh Tue Apr 7 13:25:27 GMT 2026

Reviewed & Signed 9EhQqTEH5ZQWGBNc1LEdx90D5YY

Brand Specialist - Julia Sand Tue Apr 7 18:11:40 GMT 2026

Authored & Signed EkDr5IhpNSQ5CiJXRJwIa1G2V6Q

Approval History:

Change Comment	Approved on	Status	Issued by
Updated warnings and contraindications	Tue Apr 7 2026	Approved	Abeer Prakash Singh
Updated training video	Tue Mar 17 2026	Superseded	Abeer Prakash Singh
Issue 1	Tue Mar 17 2026	Superseded	Abeer Prakash Singh
Change Comment	Approved on	Status	Issued by

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Copyright © 2026 Oasicare ApS. Alle rettigheder forbeholdes.